

# Rotate™ and Big Turn™

## Quick Reference Guide

### Operating Instructions

1. Locate both power switches on left side of control unit.
2. Turn both switches to **ON** position.



**AUTOFIRM**

### Mattress Inflation/Autofirm

Press **AUTOFIRM** button to fully inflate surface. Transfer patient to surface. **AUTOFIRM** times out after 10 minutes and resets to previous setting. To manually exit, press **AUTOFIRM** a second time.



**PATIENT  
RIGHT**



**PATIENT  
LEFT**



**FULL TURN**



**TRAINING MODE**

### Rotation Therapy

1. Choose desired turn by pressing corresponding button:
  - Right or Left Turn: maintains a turning cycle only to the selected side
  - Full Turn: provides left, center, and right rotational therapy (CLRT)
2. After selecting desired turn, press **CYCLE TIME -/+** buttons until desired time is reached. If using left or right turn as a holding cycle, press **RIGHT** or **LEFT** button until Hd shows on screen. Patient will be held on desired side.
3. After selecting therapy, pressing **TRAINING MODE** will slowly acclimate patient to therapy over a period of 60 minutes.



**PULSATE**

### Pulsation Therapy

1. Press **PULSATE** button until desired time is reached (3, 5, or 10 min).
2. Set firmness using **FIRM/SOFT** buttons based on patient comfort and clinical discretion.



**STATIC**

### Static Therapy

1. Press **STATIC** button.
2. Set firmness using **FIRM/SOFT** buttons based on patient comfort and clinical discretion.

## Other Features



FOWLER

**Fowler:** Increases airflow when head section of bed is elevated. Newer models include an Auto Fowler feature.



**Firm/Soft:** Increases (**FIRM**) or decreases airflow (**SOFT**) creating a firmer or softer surface.



LOCK OUT

**Lockout:** Locks control panel (auto after 3 min). To unlock, press and hold for 3 sec.



ALARM RESET

**Alarm Reset:** Light remains illuminated until problem is corrected.

## CPR/Quick Deflation

### 1. CPR control is on front of control unit.

Depending on model, push red button or turn dial clockwise to activate CPR. Alarm will sound.



### 2. Locate red CPR valve at head of mattress. Turn valve for quick air release.

**-OR-** Locate hose connection to control unit (all models). Disconnect for quick air release.



## Patient Transport

1. Activate Autofirm mode for 3-5 minutes to ensure mattress is completely inflated.
2. Turn both power switches off (left side of control unit). Unplug from wall.
3. Mattress will slowly deflate. As patient sinks into mattress, patient will rest on foam transport pad.

This card is not meant to replace product user manual. For complete manual, contact Sizewise.

sizewise.com 800-814-9389



**Sizewise**<sup>®</sup>  
Designed to heal