


Platinum 6000™

Quick Reference Guide



For Rev. 2 keypad, see Part No. 27650013


Sizewise®
Designed to heal

Keypad Controls



Power/Standby: Power On (Green LED); Standby (Amber LED)



Mode: Select between Therapy Modes



Autofirm: max inflate for transferring patient



Alternation Therapy



Alternation Therapy with Microclimate Management/Low Air Loss



Immersion Therapy with Microclimate Management/Low Air Loss



Immersion Therapy



Cycle Time: for Alternation Therapy only



Alarm Silence: press to reset alarms



Power Fail Warning






Low Pressure Warning



Lock/Unlock

Operating Instructions

1. Press and hold . Green LED indicates pump is on; Amber LED indicates standby mode.
2. Press  until  is selected to autofirm surface; transfer patient to surface.

Setting Therapy Mode

1. Use  to select Therapy Mode.



Alternation Therapy (no Microclimate Management/Low Air Loss)




Alternation Therapy with Microclimate Management/Low Air Loss



Immersion Therapy with Microclimate Management/Low Air Loss





Immersion Therapy (no Microclimate Management/Low Air Loss)

2. Alternation mode requires a 5, 10, or 15 min. cycle. Use  to adjust.



Patient Specific Adjustments

To customize the mattress for your patient, follow these three steps:





Patient Height: Select and use   until patient's height is shown on display.



Patient Weight: Select and use   until patient's weight is shown on display.








Comfort Settings: Select to override auto pressure settings. Use   to adjust surface's firmness (15mmHg to 50mmHg).






Changing Units of Measurement

The pump's default settings are in standard units.

Change to Metric Units:

1. Put the pump in standby by pressing .
2. Press and hold  and . Keep buttons depressed until all lights flash and you hear six beeps.
3. Cancel alarm by pressing .
4. Press  to turn on.

Change to Standard Units:

1. Put the pump in standby by pressing .
2. Press and hold  and . Keep buttons depressed until all lights flash and you hear six beeps.
3. Cancel alarm by pressing .
4. Press  to turn on.


Note: A decimal point in the Patient Weight screen indicates metric units.

CPR/Quick Deflation

1. Locate red CPR valve at head of mattress. Turn valve for quick air release.
2. Detach hoses from pump by pressing blue button and pulling connector from unit.

Patient Transport

To maintain inflation during transport:

1. Ensure CPR valve is in closed position.
2. Press  to put pump in standby mode (amber LED).
3. If mattress is to be disconnected from pump, place **Transport Cap** over hose connections.

This card is not meant to replace product user manual. For complete manual, contact Sizewise.

sizewise.com 800-814-9389

OPS-0047A-1018 27650011, Rev. 4.0

