

Looking for a Med-Surg bed with near

ZERO DOWNTIME?

Meet the new Sizewise Alliance™

Out of sight, out of mind—but still costing you money. Your hospital is not the only one with a budget-draining bed graveyard where you hide scrapped equipment that's too pricey to fix and too new to throw out. Enter Sizewise Alliance. Our new, easy-to-use, budget-friendly Med-Surg bed is available with an industry-first Shared Risk ProgramSM, ensuring maximum uptime and free loaners for three years. Pair that with our friendly 24/7/365 support staff, and you can see why Alliance is your best option now and well into the future.



Good for Hospitals, Good for Patients

- Reducing unused/underutilized beds could save your hospital 15-20 percent per year
- Sizewise Shared Risk Program ensures maximum uptime
- Equipment management solutions available for single hospitals or large health systems

Request a no-obligation bed usage review

**Call 800-814-9389 or
contact your Sizewise
Representative**

SIZEWISE

sizewise.com 800-814-9389



KEY FEATURES

- Ultra-lightweight, highly maneuverable hospital bed—no power drive needed
- Simple design is easy to clean and maintain
- Radiolucent top deck with built-in adjustable X-ray pan
- Flat deck open architecture accepts multiple mattress options—Sizewise or other brands
- At 15.25", one of the lowest Med-Surg beds on the market
- Shared Risk Program available at time of purchase

SIZEWISE

sizewise.com 800-814-9389

SPECS

BED

Low Position – Deck to Floor 15.25"
 High Position – Deck to Floor 31.5"
 Sleep Deck Width..... 35"
 Sleep Deck Length 80-88"
 Safe Working Load..... 551 lbs

SIDE RAILS

Type High Impact Plastic
 Head/Foot (4) Standard
 Long Foot Rails Optional

HEAD/FOOT BOARDS

Type High Impact Plastic
 Included..... Yes

CASTER & BRAKING SYSTEM

Lock & Steer Casters..... 6"

BED ANGLES

Backrest/Head Section..... 0-65°
 Knee Section 0-45°
 Foot Section..... 0-20°
 Trend 13°
 Reverse Trend 14°

OPTIONS & ACCESSORIES

Rest Secure System (RSS) Standard
 • Bed Exit Alarm
 • Protocol Timer
 • Weight History

Under-bed Lighting..... Standard
 Built-in X-ray Pan..... Standard
 RSS Other Languages Optional
 Long Foot Rails Optional
 Trapeze Optional
 Line Management..... Optional
 IV Pole..... Optional
 O2 Holder..... Optional
 Nurse Call Interface Cable..... Optional
 Pillow Speaker Optional
 Power Outlets (4)..... Optional
 HL7 EMR Connectivity Optional

MATTRESSES

Numerous non-powered and low air loss therapeutic surfaces available.