

Alliance™

Med-Surg Bed



Alliance is an easy-to-use, budget-friendly Med-Surg bed available with our industry-first Shared Risk ProgramSM, which ensures maximum uptime and quick repairs. It's a great solution for hospitals with budget-draining "bed graveyards." Paired with our 24/7/365 support staff, Alliance is sure to be your best option well into the future.

- *Ultra-lightweight and highly maneuverable—no Power Drive needed*
- *Simple design is easy to clean and maintain*
- *Radiolucent top deck with built-in adjustable X-ray pan*
- *At 15.25", one of the lowest Med-Surg beds available*
- *Platform auto-retracts during bed adjustment, ensuring therapeutic positioning and improved comfort*

sizewise.com 800-814-9389


Sizewise[®]
Designed to heal



- 1 — Built-in Rest Secure System™ includes bed exit alarm, protocol timer, and scale functions.
- 2 — Bed extends to accommodate mattresses up to 88”.
- 3 — Open architecture flat deck accepts multiple mattress options—Sizewise or other brands.

Key Specs

BED
 Low Position - Deck to Floor 15.25”
 High Position - Deck to Floor 31.5”
 Deck Width 35”
 Deck Length 80-88”
 Safe Working Load 551 lbs.

SIDE RAILS
 Type High Impact Plastic
 Included (4) Standard
 Long Foot Rails Optional

HEAD/FOOT BOARDS
 Type High Impact Plastic
 Included Yes

CASTER AND BRAKING SYSTEM
 Lock and Steer Casters 6”

BED ANGLES
 Backrest/Head Section 0-65°
 Knee Section 0-45°
 Foot Section 0-20°
 Trend 13°
 Reverse Trend 14°

OPTIONS AND ACCESSORIES
 Rest Secure System™ (RSS) Standard
Bed Exit Alarm, Protocol Timer, Built-in Scale
 Under-bed Lighting Standard
 Therapeutic Auto-Retract Standard
 RSS Other Languages Optional
 Trapeze Optional
 Line Management Optional
 IV Pole Optional
 O2 Holder Optional
 Nurse Call Interface Cable Optional
 Pillow Speaker Optional
 Power Outlets (4) Optional
 HL7 EMR Connectivity Optional

MATTRESS OPTIONS
 Numerous non-powered and low air loss therapeutic surfaces available.