


Quick-Ship Products

Bed and Support Surface Options



Dual-Core Mattress



Power Foundation



Dual-Core Mattress
w/Defined Perimeter Edge

Sizewise Quick-Ship options provide facilities a cost-effective way to quickly increase patient capacity. The **Dual-Core Mattress** features two layers of foam for comfort and pressure redistribution. Polycarbonate 4-Way Stretch Top Cover fabric reduces friction and shear while standing up to rigorous hospital cleaning protocols. Optional defined perimeter edge provides additional support. Fire barrier meets CFR1633.

Power Foundation features head-up/foot-up and one-touch return to flat using hand control. Head-of-bed angle from zero to 35 degrees. Three static height options available during setup. Foundation includes wipe-down cover. Flat Foundation option also available.

Delivery Time: Sizewise help deliver these products in large quantities in a very short amount of time—typically within a week. Ask your local Sizewise representative to confirm quantity and speed of delivery.

Note: the power foundation is widely used in the consumer market, and is certified to meet UL 962. This product is not certified to meet 60601 standards.

Key Specs

DUAL-CORE MATTRESS

Capacity	400 lbs.
Width	36"
Length	80"
Polycarbonate 4-Way Stretch Top Cover ..	Standard
Vinyl Laminate Bottom Cover	Standard
Fire Barrier (CFR-1633)	Standard
Edge support	Standard
Heel slope	Standard
Defined perimeter edge	Optional

POWER FOUNDATION

Capacity	850 lbs.
Width	36"
Length	80"
Static Height Options	7", 11", 15"
Certification	UL 962
Battery Back-Up	Standard
Head Angle	0-35°
Trend/Rev Trend	Not Available
Patient Hand Control	Standard
Foundation Wipe-Down Cover	Included

sizewise.com 800-814-9389

