

Stratus 5.1™ and 3.1™

Low Air Loss Alternation Therapy



Stratus 5.1™

Stratus 5.1™ and Stratus 3.1™ offer clinically effective 1-in-2 alternation therapy for the prevention and treatment of pressure injuries. Each system features 20 individual high-strength air cushions for pressure redistribution and blood flow stimulation, powered by an easy-to-use pump.

- Two therapy modes: 1-in-2 alternation therapy and static
- Variable cycle times: 10-25 min. (Stratus 5.1)
- Whisper-quiet pumps with adjustable pressure settings

Key Specs

MATTRESS

Capacity (Standard)	550 lbs.
Capacity (Bariatric)	650 lbs.
Width	35"
Length	80"
Standard and Bariatric Sizes	Available

CONTROL UNITS

Stratus 5.1 (Digital)

Liters of Air Per Minute	7
Cycle Times	10-25 min
Zoned Air Loss	Yes
Alarms	Low Pressure, Power/System Failure
Quick Inflator	Optional

Stratus 3.1 (Dial)

Liters of Air Per Minute	7
Cycle Time	10 min
Alarms	Low Pressure
Quick Inflator	Optional

size wise.com 800-814-9389


Size wise®
Designed to heal