

Horizon 13.1™

Long-Term Care Bed Frame



Horizon 13.1™ is a long-term care bed frame with an expandable/retractable deck for accommodating all patient sizes. An adjustable wall saver protects walls from bed edges, and synchronized high-low keeps the bed level for positioning quietly and smoothly.

- *Expandable/retractable deck width: 35" to 42"*
- *Adjustable wall saver to protect walls*
- *Synchronized high-low to keep bed level*
- *Easy-to-use 2-pedal caster locking system*
- *Optional trapeze kit*

Key Specs

BED

Low Position – Deck to Floor	7"
High Position – Deck to Floor	30"
Deck Width	35"-42"
Deck Length	76" or 80"
Safe Working Load	600 lbs.

SIDE RAILS

Premium Counter-Rotating Assist Device ..	Optional
Premium Rotating Assist Bar	Optional
Premium Fixed Assist Bar	Optional

HEAD/FOOT BOARDS

Shape	Square or Center Crown
Included	Yes
Color Options	3 Available Choices
<i>Solar Oak, Medium Oak, Windsor Mahogany</i>	

BED ANGLES

Trend	Standard
Reverse Trend	Standard
Chair	Standard

OPTIONS AND ACCESSORIES

Embedded Staff Control	Optional
Enhanced Embedded Staff Control	Optional
USB Power Supply	Optional
Adjustable Length (84"-88")	Optional
Trapeze Kit	Optional
Pan Deck Kit	Optional

size-wise.com 800-814-9389

