

Hercules Disposable Sheet™



Key Specs

Weight Capacity	750 lbs.
Sheet Width	42"
Sheet Length	192"
Fabric type	100% Non-Woven Poly
Recyclable	Yes

The Hercules Disposable Sheet is designed to meet the needs of today's patient and the objectives of today's healthcare providers. The sheet is made of a high performance, 100% recyclable, non-woven fabric that creates a soft and comfortable surface. Compatible with 39" Hercules™ and Air Hercules Patient Repositioner™.

- *Affordable, single-patient-use sheet eliminates the logistical challenges of special laundry*
- *Hypoallergenic and lint-free sheet provides a clean, protective layer between patient and support surface*
- *Highly breathable, silk-like fabric transfers body heat*
- *Sized to provide a perfect fit and to minimize wrinkling*
- *Compatible with 39" Hercules™ and Air Hercules Patient Repositioner™*

size-wise.com 800-814-9389

