

Alternate™ /Mighty Air™

Low Air Loss Alternation Therapy



Alternate™





Mighty Air™

Alternate™ and its bariatric counterpart, Mighty Air™, offer optimal pressure redistribution with 1-in-3 alternation therapy, targeting pressure injuries by increasing capillary blood flow to the skin. Static mode promotes proper envelopment for microclimate management. Blower provides 1,740 liters of air per minute.

- Pressure injury prevention and treatment
- 1-in-3 alternation therapy targets skin breakdown and increases capillary blood flow to the skin
- Vapor-permeable, low-friction, and anti-shear top cover
- Hill-Rom® TotalCare® and VersaCare® replacements available

sizewise.com 800-814-9389


Sizewise®
Designed to heal



Key Specs

MATTRESS

Capacity 600 lbs. or 1,000 lbs.
 Width 35", 39", or 48"
 Length 80"

CONTROL UNIT

Liters of Air Per Minute 1,740
 Zoned Low Air Loss Yes
 Alarms Power Failure and Low Pressure
 Pressure Sensor Technology Yes
 CPR Release Yes
 Unit Weight 16 lbs.

FABRIC OPTIONS

Nylon 70 Denier-Polyurethane Coated Standard
 Polycarbonate 4-Way Stretch Optional
 Non-Skid Polyurethane Bottom Cover Standard

- 1 — 1-in-3 alternation therapy provides efficient pressure redistribution
- 2 — 1,740 liters of air per minute assists in managing microclimate
- 3 — Hill-Rom® TotalCare® and VersaCare® replacements available (TotalCare shown)

size wise.com 800-814-9389



PRD-0026-0418 Hill-Rom, TotalCare, and VersaCare are trademarks of Hill-Rom. Replacement support surfaces are designed under special license.