


# Carewise 3.1™

Foam Support Surface



## Key Specs

Capacity .....	300 lbs.
Width .....	35"
Length .....	.80" or 84"
Height.....	7"

## CONSTRUCTION

Foam Layers.....	2
------------------	---

## FABRIC OPTIONS

4-Way Stretch Top Cover .....	Standard
PVC Laminate Bottom Cover .....	Standard
Fire Barrier .....	Standard

Carewise 3.1™ offers two layers of foam: softer on top with a reverse convolute-cut bottom layer to prevent bottoming out. Built-in foam perimeters aid in easy transfers and fall prevention. And the top cover is moisture-resistant, antimicrobial, and low-shear.

- *Soft foam top layer for comfort*
- *Reverse convolute-cut bottom layer for the prevention of bottoming out*
- *Built-in foam perimeters for easy transfers and fall prevention*
- *Moisture-resistant, antimicrobial, low-shear top cover*

size wise.com 800-814-9389

