


Carewise 3.1[®]

Foam Support Surface



Key Specs

Capacity	300 lb.
Width	35"
Length80" or 84"
Height.....	7"

CONSTRUCTION

Foam Layers.....	2
------------------	---

FABRIC OPTIONS

4-Way Stretch Top Cover	Standard
PVC Laminate Bottom Cover	Standard
Fire Barrier	Standard

Carewise 3.1[®] offers two layers of foam: softer on top with a reverse convolute-cut bottom layer to prevent bottoming out. Built-in foam perimeters aid in easy transfers and fall prevention. And the top cover is moisture-resistant, antimicrobial, and low-shear.

- *Soft foam top layer for comfort*
- *Reverse convolute-cut bottom layer for the prevention of bottoming out*
- *Built-in foam perimeters for easy transfers and fall prevention*
- *Moisture-resistant, antimicrobial, low-shear top cover*

size wise.com 800-814-9389

