

# SIZEWISE™: Defined Perimeter Pressure Redistribution



## Sizewise Foam Mattress

Featured on SW Lowboy® Bed

- Pressure redistribution
- Top cover – 70 denier Butyl backed nylon taffeta, 360° Zipper
- Bottom cover – Non Skid with Evacuation Handles
- 6", 1.5 lbs, 33 ILD Medical Grade Foam
- Fire Barrier
- Raised Side Bolsters with center egress

**MADE  
IN THE USA**



CONTACT SIZEWISE: 800.814.9389 | [sizewise.net](http://sizewise.net)