

Mighty Rest®

Bariatric Rehab/Fall Prevention Bed



Mighty Rest® is designed to promote a safe and comfortable environment for patients weighing up to 1,000 lbs. The bed is available in semi- and fully-electric models with head, foot, and bed height adjustability. Mighty Rest comes in four widths and is available in extra-long sizes.

- Superior comfort and functionality
- Multiple widths and lengths
- Eight-function hand control
- Includes bariatric foam mattress and head-end side rails

Key Specs

BED

Low Position – Deck to Floor	13.5"
High Position – Deck to Floor	30"
Deck Width	39", 48", 54", or 60"
Deck Length	80" or 86"
Safe Working Load	1,000 lbs.

SIDE RAILS

Head-end	Standard
----------	----------

HEAD/FOOT BOARDS

Material	High-Impact Plastic
----------	---------------------

BED ANGLES

Head Section	0-64°
Knee Section	0-54°
Foot Section	0-10°
Trend	11°
Reverse Trend	-11°

OPTIONS AND ACCESSORIES

Bariatric Foam Mattress	Standard
Hand Control (4-button)	Standard
Hand Control (8-button)	Optional
IV Pole	Optional
Trapeze	Optional

size wise.com 800-814-9389

