

Bari Drop-Arm Commode™

Bariatric Personal Care Item



Key Specs

Capacity	600/850/1,000 lbs.
Seat Width	20"/24"/28"/36"
Overall Width	24"/28"/32"/40"
Height	18-22"/19"

Bari Drop-Arm Commode™ comes in four sizes: a 20" seat width (600 lb. capacity); 28" seat width (1,000 lb. capacity); and 36" seat width (1,000 lb. capacity). Drop-Arm Commode features removable back and armrests, and is cushioned for comfort.

- Available in two sizes
- Removable back and armrests
- Cushioned for comfort
- 24" size comes with adjustable leg height

size-wise.com 800-814-9389

