

Evolution[®]

Fall Prevention Bed



Shown with optional defined perimeter foam support surface

Evolution is a go-anywhere fall prevention bed with a 9" deck low and 30" deck high for an ergonomic work environment. The bed transports patients in any high or low position. Rest Secure System™ (RSS) comes standard. Combine Evolution with a low air loss surface to provide alternation, pulsation, rotation, or immersion therapy.

- *Focus on fall prevention: 9" deck low*
- *Rest Secure System™ (RSS) includes bed exit alarm, protocol timer, and scale*
- *Transports in any high or low position*
- *Evolution Expandable™ accommodates 42" or 48" mattresses*

size wise.com 800-814-9389


Size wise[®]
Designed to heal



- 1 — The standard Evolution® has a deck low of 9", the lowest full-feature bed frame. Also available, Evolution Low Boy®, with an industry leading 7.5" low deck.
- 2 — Optional expandable deck: adjusts from 39" to 42" or 48"
- 3 — Rest Secure System™ includes bed exit alarm, protocol timer, and scale

Key Specs

BED

Low Position – Deck to Floor9"
High Position – Deck to Floor30.75"
Deck Width (Standard)35" or 39"
Deck Width (Expandable)39", 42", 48"
Deck Length82"
Safe Working Load850 lbs.

SIDE RAILS

Type Fold Down
Included 2 Standard; Foot Optional

HEAD/FOOT BOARDS

Type High-Impact Plastic
------	---------------------------

CASTER AND BRAKING SYSTEM

Caster Diameter4"
Dual Brake/Steer System4"

BED ANGLES

Head Section 0-58°
Knee Section 0-49°
Trend 15°
Reverse Trend-15°

OPTIONS AND ACCESSORIES

Battery Backup Standard
CPR Release Standard
Floor Safety Mat Standard
Footboard Lockout Controls Standard
Rest Secure System™ (RSS) Standard
<i>Bed Exit Alarm, Protocol Timer, Built-in Scale</i>	
IV Pole Optional
O2 Holder Optional
Pillow Speaker Optional
Seizure Pad Covers Optional
Trapeze Optional
6-Port Aux Outlet Optional
Color: Mist Standard
Color: Cocoa Optional

MATTRESS OPTIONS

Numerous non-powered and low air loss therapeutic surfaces available.