

# Bari Classic™

Bariatric Rehab/ICU Bed



Bari Classic™ promotes a safe and comfortable environment for patient and caregiver alike. The bed transports in any low or high position, and features Cardiac Chair, Trend, and Reverse Trend positioning. With a safe working load of 1,000 lbs. and an expandable deck, Bari Classic truly is *Designed to heal™*.

- 1,000 lbs. safe working load
- Transports in any high/low position
- Includes Rest Secure System™ with bed exit alarm, protocol timer, and scale
- Available extra long deck (86")
- Expandable deck: 39" to 48"

size-wise.com 800-814-9389

  
**Size-wise®**  
Designed to heal



1 — Expandable deck: 39" to 48"

2 — Rest Secure System™ includes bed exit alarm, protocol timer, and scale (standard on capital purchase only)

## Key Specs

### BED

Low Position – Deck to Floor .....	15.5"
High Position – Deck to Floor .....	30"
Deck Width.....	39" or 48"
Deck Length .....	80 or 86"
Safe Working Load.....	1,000 lbs.

### SIDE RAILS

Type. ....	Removable 3-Position Fold Down
Included .....	2 Standard; Foot Optional

### HEAD/FOOT BOARDS

Material .....	High-Impact Plastic
Included .....	Yes

### CASTER AND BRAKING SYSTEM

Dual Brake/Steer System .....	5"
-------------------------------	----

### BED ANGLES

Backrest/Head Section .....	0-64°
Knee Section .....	0-54°
Foot Section .....	0-10°
Trend .....	11°
Reverse Trend .....	-11°

### OPTIONS AND ACCESSORIES

Rest Secure System™ (RSS) .....	Standard
<i>Bed Exit Alarm, Protocol Timer, Built-in Scale</i>	
Battery Backup .....	Optional
IV Pole .....	Optional
O2 Holder .....	Optional
Pillow Speaker .....	Optional
Power Drive .....	Optional
Full-Frame Trapeze .....	Optional
6-Port Aux Outlet .....	Optional

### MATTRESS OPTIONS

NP12™ pressure redistribution surface included with purchase. Numerous non-powered and low air loss therapeutic surfaces available.