

Legacy/APM™ Control Unit

Alternation Therapy Air Pump



Key Specs

CONTROL UNIT

Cycle Times 10, 15, 20, 25 min.
Alarm Low Pressure
Pressure Sensor Technology Yes
Unit Weight 4.5 lbs.

MATTRESS OPTIONS

Legacy/APM™ Control Unit available with NPT3™, NP12™, and NP9™ non-powered mattresses.

The Legacy/APM Control Unit is an option for Sizewise NPT3™, NP12™, and NP9™ foam mattresses. The pump provides Static and Alternation therapy and is fully adjustable for patient comfort and healing. Simple controls include cycle times (10 to 25 minutes), soft/firm, and Autofirm max-inflate function. Control unit hangs from the footboard of most medical bed frames.

- Static and 1-in-2 Alternation therapy air pump
- 10, 15, 20, or 25-minute Alternation cycle times
- Whisper quiet with adjustable pressure settings for patient comfort
- Constant monitoring of air pressure; low pressure alarm and LED

sizewise.com 800-814-9389

