

# Behavioral Health Bed™

Fully Electric, Ligature-Resistant Bed



Featured in the  
Behavioral Health Design Guide

The first of its kind, Behavioral Health Bed™ is a fully electric bed frame designed in collaboration with behavioral health clinicians to meet the needs of dual-diagnosis patients. The bed ensures patient safety by using PowerSafe™ tamper-resistant fasteners on all electrical parts and lockout controls accessible with a universal key. An anchor point allows the bed to be secured to the floor.

- Fully electric with easy-to-use footboard lockout controls
- Electrical articulation for head, knee, and high/low systems
- Contiguous molded head/footboards and side rails
- Tamper-resistant features include encased wiring and security screws
- Designed in collaboration with behavioral health clinicians

[sizewise.com](http://sizewise.com) 800-814-9389

  
**Sizewise**®  
Designed to heal



- 1 — Tamper-resistant fasteners and security screws
- 2 — Footboard lockout controls
- 3 — Optional Rest Secure System™ with bed exit alarm, protocol timer, and scale

## Key Specs

### BED

Low Position - Deck to Floor	7.5"
Low Position - Deck to Floor (w/RSS)	9"
High Position - Deck to Floor	30"
Deck Width	35" or 39"
Deck Length	82"
Safe Working Load	850 lbs.
Bariatric Version	Available

### SIDE RAILS

Type	Fold Down
Head Section	Standard
Foot Section	Optional

### HEAD/FOOTBOARDS

Type	High-Impact Plastic
------	---------------------

### CASTER AND BRAKING SYSTEM

Caster Diameter	4"
Dual Brake/Steer System	4"

### BED ANGLES

Head Section	0-58°
Knee Section	0-49°
Foot Section	0-8°
Trend	15°
Reverse Trend	-15°

### OPTIONS AND ACCESSORIES

Battery Backup	Standard
Footboard Lockout Controls	Standard
Quick-disconnect hand control	Standard
Floor Safety Mat	Optional
Locking Storage Box Assembly	Optional
Pedal Lock Kit	Optional
Restraint Slots	Optional
Rest Secure System™ (RSS)	Optional
<i>Bed Exit Alarm, Protocol Timer, Scale</i>	
Seizure Pad Covers	Optional
Color: Mist	Standard
Color: Cocoa	Optional

### MATTRESS OPTIONS

Numerous non-powered and low air loss therapeutic surfaces available.

sizewise.com 800-814-9389

