

# Sapphire 1100EC™

Low Air Loss Therapy



## Key Specs

### MATTRESS

Capacity .....	1,000 lbs.
Widths .....	39", 42", 48", 54", or 60"
Length .....	80"

### CONTROL UNIT

Liters of Air Per Minute .....	1740
Cycle Times .....	3-20 min.
Alarms .....	Power Failure and Low Pressure

### FABRIC OPTIONS

Nylon 70 Denier Polyurethane-coated ...	Standard
4-Way Stretch .....	Optional
Non-Skid Polyurethane Bottom Cover ...	Standard

Sapphire 1100EC™ combines clinically effective low air loss with 1-in-3 alternation therapy in a complete mattress replacement system. Variable cycle times can adjust from three to 20 minutes. The mattress features high-strength air chambers, vented to provide the low air loss therapy and keep the patient dry and comfortable.

- *Two therapy modes: 1-in-3 alternation and static low air loss*
- *Variable rotation cycle times: three to 20 minutes*
- *Fowler feature boosts seat inflation when head is elevated*
- *Quilted top cover is vapor-permeable and machine-washable*

sizewise.com 800-814-9389