


Carewise 5.1™

Foam Support Surface



Key Specs

Capacity	400 lbs.
Width	35", 39", or 42"
Length	76", 80", or 84"
Height.....	7"

CONSTRUCTION

Foam Layers.....	3
------------------	---

FABRIC OPTIONS

Nylon 70 Denier PVC.....	Standard
PVC Laminate Bottom Cover	Standard
Fire Barrier	Standard

Carewise 5.1™ features a diamond-cut top foam layer for pressure redistribution, shear reduction, and air circulation, keeping the surface cool for optimal patient comfort. The high-resiliency heel area is softer than the rest of the surface to protect vulnerable heels, and the full-length base prevents bottoming out.

- *Diamond-cut top layer to keep patient cool*
- *High-resiliency heel area is softer than rest of surface*
- *Full-length base to prevent bottoming out*
- *Low-shear, antimicrobial, fluid-resistant top cover*
- *Made-to-order item*

size wise.com 800-814-9389