

Advantage Wheelchair™

Bariatric Mobility Item



Key Specs

Capacity 850 lbs.
Seat Width 24", 26", 28", 30", 32", and 34"
Seat Depth 20"

Advantage Wheelchair™ is constructed with a heavy-duty steel frame, four-point telescoping seat, and signature double x-brace support. Swing-away, removable footrests come standard, and the wheelchair folds for easy transport. Available in multiple seat widths.

- *Heavy-duty steel frame*
- *Four-point telescoping seat*
- *Swing-away, removable footrests*
- *Folds for easy transport*
- *Multiple seat widths available*

size-wise.com 800-814-9389


Size-wise®
Designed to heal