

# Victory 1300™

Bariatric Foam Support Surface



Victory 1300™ foam support surface provides pressure redistribution and shear reduction for optimal patient comfort using four layers of impression load deflection (ILD) foam. High-resiliency foot end is softer than the rest of the surface to protect vulnerable heels. Concealed zipper prevents skin tears and protects inner foam from fluid ingress. Optional Mattress Asset Protection (MAP) Liner™ safeguards internal components and reduces shear.

- *Low-shear, antimicrobial, fluid-resistant top cover*
- *Optional RF-welded seams for max infection control*
- *Made-to-order item*

## Key Specs

Capacity	800 lbs.
Width	48"
Length	80" or 86"
Height	7"

## CONSTRUCTION

Foam Layers	4
MAP™ Liner	Standard
Fire Barrier	Standard
Defined Perimeter	Optional

## FABRIC OPTIONS

Envy® 4-Way Stretch Top Cover	Standard
Envy® Non-Skid Polyurethane	Standard
Nylon 70 Denier PVC	Optional
RF-Welded Seams	Optional
Polycarbonate 4-Way Stretch Top Cover	Optional

sizeWise.com 800-814-9389

