

Horizon 9.1™

Long-Term Care Bed Frame



Horizon 9.1™ is a long-term care bed frame with synchronized high-low for positioning quietly and smoothly, as well as sealed ball bearings in the joints to provide continuous, silent operation. An adjustable wall saver protects walls from bed edges.

- Synchronized high-low to keep bed level
- Adjustable wall saver to protect walls
- Easy-to-use 2-pedal caster locking system
- Optional trapeze kit
- Optional deck extension kit

Key Specs

BED

Low Position – Deck to Floor	7"
High Position – Deck to Floor	30"
Deck Width	35"
Deck Length	76" or 80"
Safe Working Load	600 lbs.

SIDE RAILS

Premium Counter-Rotating Assist Device ..	Optional
Premium Rotating Assist Bar	Optional
Premium Fixed Assist Bar	Optional

HEAD/FOOT BOARDS

Shape	Square or Center Crown
Included	Yes
Color Options	3 Available Choices
<i>Solar Oak, Medium Oak, Windsor Mahogany</i>	

BED ANGLES

Trend	Optional
Reverse Trend	Optional
Chair	Optional

OPTIONS AND ACCESSORIES

Embedded Staff Control	Optional
Adjustable Width (35"-42")	Optional
Adjustable Length (84"-88")	Optional
Trapeze Kit	Optional
Pan Deck Kit	Optional

size-wise.com 800-814-9389


Size-wise®
Designed to heal