

Geri Air™

Therapeutic Chair Overlay



Key Specs

MATTRESS

Capacity	250 lbs.
Width	22"
Length	72"

CONTROL UNIT

Cycle Time	10 min.
------------------	---------

FABRIC OPTIONS

4-Way Stretch	Standard
Polycarbonate 4-Way Stretch	Optional
Non-Skid PVC Bottom Cover	Standard

Geri Air™ provides 2 in 1 alternation therapy for geriatric patients using a bedside chair. The Geri Air pad is covered with an all-way stretch, quilted nylon comforter designed to prevent friction and shear. 8 liters of air per minute.

- 2 in 1 alternation therapy
- Nylon comforter to prevent friction and shear
- Quiet control unit

size wise.com 800-814-9389