


# Crib Mattress

Foam Support Surface



## Key Specs

### MATTRESS

Capacity .....	100 lbs.
Width .....	Custom
Length .....	Custom

### FABRIC OPTIONS

Polycarbonate 4-Way Stretch .....	Standard
4-Way Stretch (in Royal or Green) .....	Optional
Non-Skid Polyurethane Bottom Cover ...	Standard

Our Crib Mattress is latex-free and radiolucent for X-rays. The First Point of Contact™ top cover is free from chemicals of concern and is vapor-permeable, low-friction, and anti-shear. Crib mattresses are available with foam or tri-laminate foam with either a memory foam or gel top layer. Available in standard and custom sizes.

- *Latex-free and radiolucent*
- *First Point of Contact™ top cover free of chemicals of concern, vapor-permeable, low-friction, and anti-shear*
- *Halogen-free fire barrier*
- *Available with foam or tri-laminate foam*
- *Available with either memory foam or gel top layer*
- *Standard and custom sizes to fit any crib or bed*

size-wise.com 800-814-9389