

Envy Line® G Series™

Non-Powered Mattress




NP12™ shown with optional green top cover and alternating pump.

Envy Line® G Series™ mattresses feature F3 Free™ construction, reducing harmful chemicals in the fabric, foam, and fire barrier. They conform readily to any surface, maximizing patient immersion while reducing shear. They are available with 9 or 12 independently conforming chambers for effective pressure redistribution and pressure injury prevention. G Series mattresses are recyclable at the end of their usable lives.

- F3 Free construction reduces harmful chemicals in the fabric, foam, and fire barrier
- Three unique zones deliver proper support and pressure redistribution
- Available as non-powered or with 1-in-2 alternation therapy
- 9 or 12 independently conforming chambers (NP9™ and NP12™)

sizewise.com 800-814-9389


Sizewise®
Designed to heal



Key Specs

MATTRESS

Capacity500 lbs. / 1,000 lbs.
 Width 35", 39", 42", or 48"
 Length80"-86"
 Custom Sizes Available

FABRIC OPTIONS

4-Way Stretch (Royal) Top Cover Standard
 4-Way Stretch (Green) Top Cover Optional
 Polycarbonate 4-Way Stretch Top Cover Optional
 Non-Skid Polyurethane Bottom Cover Standard
 Halogen-Free Fire Barrier Standard
 MAP Liner™ Optional



- 1 — Pressure Mapping: 12-chamber mattress
- 2 — Hill-Rom TotalCare® and VersaCare® replacements available
- 3 — Available with 1-in-2 alternation therapy



A PRACTICE GREENHEALTH PROGRAM

Sizewise proudly supports Practice Greenhealth's Healthier Hospitals Safer Chemicals Challenge for Healthy Interiors. Our Envy Line products meet or exceed Healthier Hospitals goals of reducing patient and caregiver exposure to chemicals of concern.

sizewise.com 800-814-9389

