

# Envy® G Series™

Non-Powered Pressure Redistribution Mattress




NP12™ shown with  
alternating pump.

Envy® G Series™ mattresses feature F3 Free™ construction, reducing harmful chemicals in the fabric, foam, and fire barrier. They conform readily to any surface, maximizing patient immersion while reducing shear. They are available with 9 or 12 independently conforming chambers for effective pressure redistribution and pressure injury prevention. G Series mattresses are recyclable at the end of their usable lives.

- F3 Free construction reduces harmful chemicals in the fabric, foam, and fire barrier
- Three unique zones deliver proper support and pressure redistribution
- Available as non-powered or with 1-in-2 alternation therapy
- 9 or 12 independently conforming chambers (NP9™ and NP12™)

size wise.com 800-814-9389

  
**Size wise**®  
Designed to heal



## Key Specs

### MATTRESS

Capacity .....500 lbs. / 1,000 lbs.  
 Width ..... 35", 39", 42", or 48"  
 Length .....80"-86"  
 Custom Sizes ..... Available

### FABRIC OPTIONS

Envy® 4-Way Stretch Top Cover .....Standard  
 Extreme Polycarb 4-Way Stretch Top ..... Optional  
 Non-Skid Polyurethane Bottom Cover .....Standard  
 Halogen-Free Fire Barrier.....Standard  
 MAP Liner™ ..... Optional

### FABRIC CHARACTERISTICS



4-Way  
Stretch



Durability



Green



1 — Pressure Mapping: 12-chamber mattress

2 — Hillrom™ TotalCare® and VersaCare® replacements available

3 — Available with 1-in-2 alternation therapy



A PRACTICE GREENHEALTH PROGRAM

Sizewise proudly supports Practice Greenhealth's Healthier Hospitals Safer Chemicals Challenge for Healthy Interiors. Our Envy products meet or exceed Healthier Hospitals goals of reducing patient and caregiver exposure to chemicals of concern.

sizewise.com 800-814-9389

