

# Horizon 3.1™

Long-Term Care Bed Frame



Horizon 3.1™ is a long-term care bed frame with synchronized high-low for positioning quietly and smoothly. An adjustable wall saver protects walls from bed edges.

- *Synchronized high-low to keep bed level*
- *Adjustable wall saver to protect walls*
- *Optional deck extension kit*
- *Optional trapeze kit*

## Key Specs

### BED

Low Position – Deck to Floor	9.95"
High Position – Deck to Floor	24.4"
Deck Width	35"
Deck Length	80-84"
Safe Working Load	450 lbs.

### SIDE RAILS

Premium Counter-Rotating Assist Device	Optional
Premium Basic Assist Bar	Optional
Premium Fixed Assist Bar	Optional

### HEAD/FOOT BOARDS

Shape	Square or Center Crown
Included	Optional
Color Options	5 Available Choices

*Solar Oak, Medium Oak, Wild Cherry,  
Montana Walnut, Windsor Mahogany*

### OPTIONS AND ACCESSORIES

USB Power Supply	Optional
Trapeze Kit	Optional
Pan Deck Kit	Optional

sizewise.com 800-814-9389