

# Bari Walker™

Bariatric Mobility Item



## Key Specs

Capacity ..... 750 lbs.  
Width Between Supports ..... 27"  
Overall Width ..... 35.5"  
Height ..... Adjustable from 31" to 40"

At 27" between supports, Bari Walker™ is wider than the average walker to fit bariatric patients. The walker is height-adjustable from 31" to 40" and features a weight capacity of 750 lbs. Folds to 5" width for easy transport.

- 27" width between supports
- Height-adjustable from 31" to 40"
- Weight capacity: 750 lbs.
- Folds to 5" width for easy transport

size-wise.com 800-814-9389

