





INFECTION CONTROL: ARE YOU AT RISK?

SizeWise Solutions: Mattress Asset Protection™ Liner

The SizeWise Mattress Asset Protection™ liner is a safe, simple and secure method of protecting medical mattresses from contamination risks in order to safeguard against patient infection.

MAP Liner is PVC Free



- The SizeWise MAP Liner is an inner liner that is below the mattress top cover, but above the foam construction.
 - Its unique feel is soft and light like silk, but is strong and durable against any seepage that may occur from mattress wear or tear.
 - An additional benefit to the patient is shear reduction - the two surfaces slide easily with one another, and work together to reduce shear.
- The New SizeWise Mattress Asset Protection™ or MAP Liner is a polyurethane liner designed specifically to address concerns from the FDA communication on April 19, 2013, stating “Damaged or Worn Covers for Medical Mattresses Pose Risk of Contamination and Patient Infection.”

