

Envy® O Series™

Pressure Redistribution OR Table Pad





Envy O Series features F3 Free™ construction that reduces harmful chemicals in the fabric, foam, and fire barrier. A four-way stretch top cover allows immersion into the pad for optimal pressure redistribution in order to prevent pressure injuries. The universal design is radiolucent and fits most Operating Room tables.

- F3 Free construction reduces harmful chemicals in the fabric, foam, and fire barrier
- Optimal immersion for pressure redistribution and PI prevention
- Five-piece universal design fits most Operating Room tables
- Custom sizes available

sizewise.com 800-814-9389


Sizewise®
Designed to heal



1 — Polymer gel topper contours to patient body while providing low-interface pressure

2 — EcoEase™ pieces are easily replaceable by section

Key Specs

SUPPORT SURFACE

Capacity	.500 lb.
Head	20"W x 10.5"L
Torso	20"W x 42.375"L
Leg	20"W x 18.5"L
Arm	6"W x 26"L
Custom Sizes	Available

FABRIC OPTIONS

Extreme Polycarb 4-Way Stretch Top	Standard
Envy® Non-Skid Bottom Cover	Standard
Halogen-Free Fire Barrier	Standard
MAP Liner™	Optional

FABRIC CHARACTERISTICS



4-Way
Stretch



Chemical
Resistant



Green

sizewise.com 800-814-9389



SizeWise®
Designed to heal