

# Envy® O Series™

Pressure Redistribution OR Table Pad




Envy O Series features F3 Free™ construction that reduces harmful chemicals in the fabric, foam, and fire barrier. A four-way stretch top cover allows immersion into the pad for optimal pressure redistribution in order to prevent pressure injuries. The universal design is radiolucent and fits most Operating Room tables.

- *F3 Free construction reduces harmful chemicals in the fabric, foam, and fire barrier*
- *Optimal immersion for pressure redistribution and PI prevention*
- *Five-piece universal design fits most Operating Room tables*
- *Custom sizes available*

sizewise.com 800-814-9389

  
**Sizewise®**  
Designed to heal



- 1 — Polymer gel topper contours to patient body while providing low-interface pressure
- 2 — EcoEase™ pieces are easily replaceable by section

## Key Specs

### BED

Capacity .....	500 lbs.
Head .....	20"W x 10.5"L
Torso .....	20"W x 42.375"L
Leg .....	20"W x 18.5"L
Arm .....	6"W x 26"L
Custom Sizes .....	Available

### FABRIC OPTIONS

Extreme Polycarb 4-Way Stretch Top .....	Standard
Envy® Non-Skid Bottom Cover .....	Standard
Halogen-Free Fire Barrier .....	Standard
MAP Liner™ .....	Optional

### FABRIC CHARACTERISTICS



4-Way  
Stretch



Chemical  
Resistant



Green



A PRACTICE GREENHEALTH PROGRAM

Sizewise proudly supports Practice Greenhealth's Healthier Hospitals Safer Chemicals Challenge for Healthy Interiors. Our Envy products meet or exceed Healthier Hospitals goals of reducing patient and caregiver exposure to chemicals of concern.

sizewise.com 800-814-9389

