

# Quick-Ship Recliner

Comfort and Mobility Item



Shown in Micron (left) and Creme fabric.\*

The Sizewise Quick-Ship Recliner is designed for both comfort and safety—allowing facilities a cost-effective way to quickly scale up patient capacity across inpatient or outpatient settings. The recliner’s upholstery and armrests are durable and easy to clean, supporting the most rigorous infection control and prevention protocols.

Quick-Ship Recliner is a great solution for acute patient rooms, emergency department overflow, outpatient vaccination centers, and other high-use areas.

- 650 lbs. weight capacity
- Durable polyurethane armrests
- Infinite position locking backrest, including full flat position
- Easy-to-use lock and steer casters
- Continuous footrest provides added comfort
- Foley bag holder

\*Note: fabric color may vary. Please contact your Sizewise rep for availability.

## Key Specs

Capacity	650 lbs.
Seat Width	24"
Seat Depth	21"
Seat Height	21.5"
Back Height	26"
Arm Height	28"

## OPTIONS

Locking Casters	Standard
Foley Bag Holder	Standard
Polyurethane Armrests	Standard
Removable Back	Standard

## DELIVERY TIME

Sizewise can help deliver these products in large quantities in a very short amount of time—typically within a week. Ask your local Sizewise representative for a Quick-Ship Quote to confirm quantity, fabric color, and speed of delivery.

[sizewise.com](https://www.sizewise.com) 800-814-9389

