

Bari Rehab Platform 3™

Bariatric Med-Surg/ICU Bed



Shown with optional pressure redistribution mattress

Bari Rehab Platform 3™ (BRP3) debuts the Power Drive 360™ for easy navigation through crowded hospital corridors—forward, backward, and side-to-side, plus the ability to turn on a dime. Available with our industry-first Shared Risk ProgramSM for maximum uptime.

- *Optional Power Drive 360™ for easy navigation*
- *Open architecture flat deck accepts most mattresses—Sizewise or other brands*
- *Fully compliant with latest IEC safety standards (60601-2-52)*
- *HL7 EMR connectivity*
- *Shared Risk Program available at time of purchase*

sizewise.com 800-814-9389


Sizewise[®]
Designed to heal



- 1 — Optional Power Drive 360™ moves forward, backward, and side-to-side with ease
- 2 — Rest Secure System 3™ includes bed exit alarm, protocol timer, and scale
- 3 — Open architecture flat deck accepts most mattresses—Sizewise or other brands
- 4 — Available with Hercules Patient Repositioner™

Key Specs

BED

Low Position – Deck to Floor	15.35"
High Position – Deck to Floor	30"
Deck Width	39-48"
Deck Length	80 or 86"
Safe Working Load	1,000 lbs.
Expandable Deck	Yes

SIDE RAILS

Type	Removable 2-Position Fold Down
Included	2 Standard; Foot Optional

HEAD/FOOT BOARDS

Type	Kydex®
Accessory Bar	Standard
Integrated Info Holder	Standard

CASTER AND BRAKING SYSTEM

Dual Brake/Steer System	5"
-------------------------	----

BED ANGLES

Backrest	0-64°
Knee Section	0-54°
Trend	15°
Reverse Trend	-15°

OPTIONS AND ACCESSORIES

Rest Secure System 3™ (RSS)	Standard
Bed Exit Alarm, Protocol Timer, Scale	
RSS Other Languages	Optional
Battery Backup	Optional
HL7 EMR Connectivity	Optional
IV Pole	Optional
Line Management	Optional
Nurse Call Interface	Standard
O2 Holder	Optional
Pillow Speaker	Optional
Trapeze	Optional
6-Port Aux Outlet	Optional

MATTRESS OPTIONS

Numerous non-powered and low air loss therapeutic surfaces available.

Hercules Patient Repositioner	Optional
Air Hercules Patient Repositioner	Optional

sizewise.com 800-814-9389

