

# Behavioral Health Bed™ – Platform

Ligature-Resistant Bed



Available in Mist (shown)  
or Cocoa.

Behavioral Health Bed™ – Platform provides optimal patient and caregiver safety. Made from high-density polyethylene and built with tamper-resistant fasteners, the Platform is durable and easy to clean. Anchor points fasten the bed to the floor, ensuring the bed is stable and secure. Available in mist or cocoa.

- *High-density polyethylene: durable and easy to clean*
- *Tamper-resistant fasteners*
- *Anchor points allow bed to be secured to floor*
- *Safe working load of 1,000 lbs.*
- *Full-enclosure base and low-profile head/footboards available*

sizewise.com 800-814-9389

  
**Sizewise**®  
Designed to heal



## Key Specs

**BED**

Sleep Surface Width .....	35"
Sleep Surface Length .....	82"
Platform Height .....	14"
Safe Working Load .....	1,000 lbs.
Product Weight (open base) .....	250 lbs.


*Note: Low-Profile Head/Footboards reduce bed weight by approximately 40 lbs.*



**OPTIONS**

Anchor Points, Floor .....	Standard
Anchor Points, Wall* .....	Optional
Colors .....	Mist or Cocoa
Custom Head/Footboard .....	Optional
Full-Enclosure Base .....	Optional
Low-Profile Head/Footboard .....	Optional
Restraint Slots .....	Optional

\*Custom anchor points available upon request and approval by Sizewise.



**MATTRESSES**  
 Sizewise Behavioral Health Mattress™ available with or without defined perimeter edge.

- 1 — Two color options: mist or cocoa (shown)
- 2 — Four anchor points allow bed to be secured to floor
- 3 — Available Behavioral Health Mattress™ designed with safety in mind
- 4 — Head/Footboard available in Standard and Low profiles

sizewise.com 800-814-9389

