

Alliance™

Med-Surg Bed




Alliance is an easy-to-use, budget-friendly Med-Surg bed available with our industry-first Shared Risk ProgramSM, which ensures maximum uptime and quick repairs. Lightweight and highly maneuverable, its simple design is easy to clean and maintain. And the platform auto-retracts during bed adjustment, ensuring therapeutic positioning and improved patient comfort.

- *Ultra-lightweight and highly maneuverable—no Power Drive needed*
- *Simple design is easy to clean and maintain*
- *Radiolucent top deck with built-in adjustable X-ray pan*
- *Available with our Shared Risk ProgramSM*
- *Therapeutic Auto-Retract ensures patient comfort during bed adjustment*

sizewise.com 800-814-9389


Sizewise[®]
Designed to heal



- 1 — Built-in Rest Secure System™ includes bed exit alarm, protocol timer, and scale functions
- 2 — Bed extends to accommodate mattresses up to 88"
- 3 — Open architecture flat deck accepts multiple mattress options—Sizewise or other brands

Key Specs

BED

Low Position - Deck to Floor	17.375"
High Position - Deck to Floor	31.5"
Deck Width	35"
Deck Length	.80-88"
Safe Working Load	550 lbs.

SIDE RAILS

Type	High Impact Plastic
Included (4)	Standard
Long Foot Rails	Optional

HEAD/FOOT BOARDS

Type	High Impact Plastic
Included	Yes

CASTER AND BRAKING SYSTEM

Lock and Steer Casters	6"
------------------------	----

BED ANGLES

Backrest/Head Section	0-65°
Knee Section	0-45°
Foot Section	0-20°
Trend	12°
Reverse Trend	-12°

OPTIONS AND ACCESSORIES

Rest Secure System™ (RSS)	Standard
<i>Bed Exit Alarm, Protocol Timer, Built-in Scale</i>	
Under-bed Lighting	Standard
Therapeutic Auto-Retract	Standard
RSS Other Languages	Optional
Trapeze	Optional
Line Management	Optional
IV Pole	Optional
O2 Holder	Optional
Nurse Call Interface Cable	Optional
Pillow Speaker	Optional
Power Outlets (4)	Optional
HL7 EMR Connectivity	Optional

MATTRESS OPTIONS

Numerous non-powered and low air loss therapeutic surfaces available.

sizewise.com 800-814-9389

