

# Seizure Pads

Safety Accessory



## Key Specs

Size ..... Varies  
Fabric ..... Envy™ 4-Way Stretch

Sizewise Seizure Pads effectively cover bed frame side rails and head/footboards, protecting patients from potential injury. The pads are made with closed-cell foam and Envy™ 4-Way Stretch fabric, ensuring a tight fit. Seams are RF-welded for durability and infection prevention. Seizure Pads are available for our Alliance™, Behavioral Health Bed™, and Evolution™ bed frames.

- Shields patient from potential injury during a seizure
- RF-welded seams for easy cleaning and infection prevention
- Available for Alliance™, Behavioral Health Bed™, and Evolution™ bed frames

[sizewise.com](https://sizewise.com) 800-814-9389

