

Advantage™

LTC and Homecare Bed




Shown with optional expandable deck and NP12 support surface

Advantage Bed promotes a safe, comfortable healing environment in Long-Term Care and home settings. Constructed with durable, heavy-duty steel, this fully electric bed has a safe working load of 850 lbs. Comes in two standard widths (39" and 48"), while Advantage Expandable™ transforms to 48" wide.

- *Fully electric bed includes Comfort Chair position*
- *Multi-function hand control and footboard lockout*
- *850 lbs. capacity*
- *Advantage Expandable adjusts to 48" wide*

size-wise.com 800-814-9389


Size-wise®
Designed to heal



- 1 — Advantage Expandable™ adjusts from 39" to 48"
- 2 — Fully electric bed includes Comfort Chair position
- 3 — Footboard lockout controls

Key Specs

BED

Low Position – Deck to Floor	14"
High Position – Deck to Floor	27"
Deck Width	39", 48", or 39-48"
Deck Length	80"
Safe Working Load	850 lbs.

SIDE RAILS

Head-end	Standard
Foot-end	Optional

HEAD/FOOT BOARDS

Material	High-Impact Plastic
----------	---------------------

BED ANGLES

Head Section	0-64°
Knee Section	0-54°
Foot Section	0-10°
Trend	11°
Reverse Trend	-11°

OPTIONS AND ACCESSORIES

Footboard Controls	Standard
Hand Control (8-button)	Standard
Battery Backup	Optional
Cord Wrap	Optional
Expandable Deck (48"W)	Optional
Foley Bag Hooks	Optional
Foot-end Side Rails	Optional
IV Pole	Optional
Manual CPR Pull	Optional
O2 Holder	Optional
Trapeze	Optional

MATTRESS OPTIONS
 Numerous non-powered and low air loss therapeutic surfaces available.